

# User Manual

## Trinocular Compound Microscope

Model M837L1



[MicroscopeNet.com](http://MicroscopeNet.com)

# Table of Contents

|                                |   |
|--------------------------------|---|
| i. Caution.....                | 1 |
| ii. Care and Maintenance.....  | 2 |
| 1. Component Illustration..... | 3 |
| 2. Installation .....          | 4 |
| 3. Operation.....              | 5 |
| 4. Specifications.....         | 7 |
| 5. Troubleshooting Guide.....  | 8 |


## i. Caution

1. Find the "UP" sign and place the Styrofoam container on your table or bench so that the arrow is pointing upward. Open the shipping carton carefully to prevent any accessory items (i.e. objectives or eyepieces) from dropping and being damaged.
2. Do not discard the molded Styrofoam container. The container should be retained should the microscope ever requires reshipment.
3. Keep the instrument out of direct sunlight, high temperature or humidity, and dusty environments. Ensure that the microscope is located on a smooth, level and firm surface.
4. If any specimen solutions or other liquids splash onto the stage, objective or any other component, disconnect the power cord immediately and wipe up the spillage. Otherwise, the instrument may be damaged.
5. All electrical connectors (power cord) should be inserted into an electrical surge suppressor to prevent damage due to voltage fluctuations.
6. Confirm that the input voltage indicated on your microscope corresponds to your line voltage. The use of a different input voltage other than that as indicated will cause severe damage to the microscope.

## ii. Care and Maintenance

1. Do not attempt to disassemble any component including eyepieces, objectives or focusing assembly.
2. Keep the instrument clean; remove dirt and debris regularly. Accumulated dirt on metal surfaces should be cleaned with a damp cloth. More persistent dirt should be removed using a mild soap solution. **Do not use organic solvents for cleansing.**
3. The outer surface of the optics should be inspected and cleaned periodically using an air stream from an air bulb. If dirt remains on the optical surface, use a soft cloth or cotton swab dampened with a lens cleaning solution (available at camera stores). All optical lenses should be swabbed using a circular motion. A small amount of absorbent cotton wound on the end of a tapered stick makes a useful tool for cleaning recessed optical surfaces. Avoid using an excessive amount of solvents as this may cause problems with optical coatings or cemented optics or the flowing solvent may pick up grease making cleaning more difficult. Oil immersion objectives should be cleaned immediately after use by removing the oil with lens tissue or a clean, soft cloth.
4. Store the instrument in a cool, dry environment. Cover the microscope with the dust cover when not in use.

# 1. Components Illustration



## 2. Installation

### 2.1 Installation of the trinocular viewing head

- 1) Loosen the thumb screw on the top of the body and remove the plastic cover on the top.
- 2) Remove the cap on the dovetail of the Trinocular viewing head.
- 3) Insert the Trinocular viewing head into the head of the body; ensure that the dovetail is completely seated into the body; tighten the thumb screw.

**Caution:** Do not release the head from your hand grip until you are sure the head is installed securely.

### 2.2 Installation of the eyepieces

- 1) Remove the protective caps from the eyepiece tubes.
- 2) Insert the eyepieces into the eyepiece tubes.

### 2.3 Installation of the objectives

- 1) Turn the coarse focus knob until the mechanical stage is at its lowest position.
- 2) Install the 4X objective into the nosepiece. Then in a clock-wise direction, rotate the nosepiece and install each succeeding higher magnification objective.

**Note:** When changing the objective magnification, rotate the objective nosepiece until you hear a “click” sound. This ensures the objective is centered in the optical light path.

### 2.4 Installing (or changing) the batteries

- 1) Turn the power off and disconnect AC adapter.
- 2) Turn over the microscope on its side; find the battery housing at the bottom.
- 3) Open the cover of the battery housing by loosening the thumb screw. Take out the dead batteries and insert the new batteries.
- 4) Put the cover on.

**Caution:** Before you turn over the microscope, be sure to take the eyepieces off and be certain that the head is securely locked by the thumb screw.



### 2.5 Replacing the fuse

- 1) Turn off the power and disconnect the AC adapter.
- 2) Turn over the microscope on its side; find the fuse holder at the bottom of the base.
- 3) Turn the fuse holder counter-clockwise to take it off, insert new fuse, and then turn it on clockwise.

**Caution:** Before you turn over the microscope, be sure to take the eyepieces off and be certain that the head is securely locked by the thumb screw.

### 3. Operation

#### 3.1 Adjusting illumination

- 1) Connect the AC adapter to the power outlet and the microscope (or install 3 AA batteries into the battery housing).
- 2) Turn on the main power switch.
- 3) Rotate the variable intensity dial to increase or decrease the brightness.


#### 3.2 Placing specimen

- 1) Place the slide on the mechanical stage. Use the slide holder to gently secure the slide.
- 2) Turn the X and Y translational knobs to position the specimen for viewing.

**Caution:** Be sure not to allow an objective to touch a specimen slide when changing objectives.

#### 3.3 Adjusting interpupillary distance

While observing with both eyes, hold the left and right eyepiece tubes then slowly slide the tubes in and out. The interpupillary distance is correct when the left and right fields of view converge completely into one image.



#### 3.4 Adjusting eyepiece diopter

- 1) Rotate the 10x objective into position.
- 2) Rotate the diopter ring on the right eyepiece tube until its numerical value is the same as your interpupillary distance, for example, 65 in the right figure.
- 3) Close your left eye and bring the specimen into focus following the focusing procedures in 3.5.
- 4) Close your right eye and bring the same specimen into clear sharp focus by adjusting only the diopter ring on left eyepiece tube. Don't use focus knobs at this step.
- 5) Since both sides are adjustable, you may also do the above in the opposite way, in other words, left eye first and right eye second.



#### 3.5 Focusing

- 1) With the 10x objective in position, raise the mechanical stage using the coarse focus knob until the specimen is close to the objective.
- 2) Turn the coarse focus knob until the specimen is in focus.
- 3) Use the fine focus knob to obtain a sharp image.
- 4) You may now switch to another magnification objective.

**Tips:** To prevent your specimen slide from making contact with an objective, raise the stage to its highest position without contacting the 100x objective, then tighten the focus stop lever.

#### 3.6 Adjusting condenser

- 1) Turn the condenser focus knob to raise or lower the condenser.

- 2) The condenser is raised when using high magnification objectives and lowered when using low magnification objectives.

**Note:**

- The centering of the condenser and the light axis of the objective are factory adjusted. Do not attempt re-adjustment.
- The highest position of the condenser has been factory adjusted. Do not attempt re-adjustment.

### 3.7 Adjusting aperture iris diaphragm


Move the aperture iris diaphragm Lever left or right to adjust the aperture size.

**Note:** The iris diaphragm is designed to adjust the aperture size, not to adjust brightness. Generally, opening the diaphragm to 70-80% of the N.A. value of the respective objective will provide an image of acceptable quality. If you want to observe the image of the iris diaphragm, remove one eyepiece and look through the tube. You will see a dark circle encroaching on the bottom of the tube.



### 3.8 Adjusting focus tension

The focus tension has been pre-set at the factory. If the mechanical stage drops by itself, rotate the tension adjustment ring situated between the coarse focus knob and microscope body on the power switch side until the tension is maintained.



### 3.9 Photo/video observing, capturing and recording

- 1) Remove the plastic photo tube port cover on the top of the viewing head.
- 2) Thread the photo tube on to the trinocular viewing head.
- 3) Turn on the computer and launch the observing software to examine.
- 4) If necessary, adjust the height by loosening the 2 set screws on the photo tube and turn the upper part in order to make the camera parfocal with the eyepieces.
- 5) You also can capture images or record live videos through the software.



**Note:** Camera and digital camera sold separately.

## 4. Specifications

### General

|                       |   |
|-----------------------|---|
| Objective Tube Length | 160mm   |
| Viewing Head          | Trinocular, inclined 45°, swiveling 360°, interpupillary distance 55-75mm, adjustable diopter |
| Eyepieces             | WF10X, P16X   |
| Nosepiece             | Quadruple   |
| Objectives            | Achromatic 4X, 10X, 40X(spring), 100X(spring, oil)  |
| Focus system          | Coaxial coarse and fine focusing, minimum fine focusing adjustment at 0.002mm, range 28mm     |
| Condenser             | Abbe, NA=1.25   |
| Stage                 | Double layer mechanical stage, area 140mmx140mm, movement range 75mm X 50mm                   |
| Illumination          | LED light   |
| Power                 | AC adapter: 100-240VAC, 50/60HZ, output DC7.5V<br>DC Battery: 3 AA batteries                  |

### Eyepieces

| Designation | Magnification | Field of View | Focal Length |
|-------------|---------------|---------------|--------------|
| Wide Field  | 10X           | 18mm          | 24.94mm      |
| Plan Field  | 16X           | 11mm          | 15.58mm      |

### Objectives

| Magnification | Numerical Aperture | Thickness of slide cover | Focus   | Working Distance | Working Mode |
|---------------|--------------------|--------------------------|---------|------------------|--------------|
| 4X            | 0.10               | 0.17mm                   | 31.04mm | 37.5mm           | Dry          |
| 10X           | 0.25               | 0.17mm                   | 17.13mm | 7.316mm          | Dry          |
| 40X(S)        | 0.65               | 0.17mm                   | 4.65mm  | 0.632mm          | Dry          |
| 100X(S, Oil)  | 1.25               | 0.17mm                   | 2.906mm | 0.198mm          | Oil          |

### Magnifications

| Eyepiece      | 10X | 16X | 10X  | 16X  | 10X  | 16X  | 10X   | 16X   |
|---------------|-----|-----|------|------|------|------|-------|-------|
| Objective     | 4X  |     | 10X  |      | 40X  |      | 100X  |       |
| Magnification | 40X | 64X | 100X | 160X | 400X | 640X | 1000X | 1600X |

## 5. Troubleshooting Guide

### OPTICAL PROBLEMS

| <b>Problem</b>  | <b>Cause</b>  | <b>Solution</b>  |
|---|---|--|
| Darkness at the periphery or uneven brightness in the field of view | Revolving nosepiece not in click stop position                                      | Revolve the nosepiece to click-stop position by swinging the objective correctly into the optical path |
|   | The light source of the bulb is not at the center                                   | Adjust the position of the bulb  |
| Dirt or dust on the view  | Dirt or dust on the lens eyepiece, condenser, objective, collector lens or specimen | Clean the lens with a camera cleaning kit  |
| Poor image quality  | No slide cover attached to the slide  | Attach a 0.17mm slide cover  |
|   | Slide cover is too thick or thin  | Use a slide cover of the appropriate thickness (0.17mm)  |
|   | Slide may be upside down  | Turn slide over so the cover-glass faces up  |
|   | Immersion oil is on a dry objective (especially the 40x)                            | Check the objectives, clean if necessary   |
|   | No immersion oil used with 100x objective   | Use immersion oil  |
|   | Air bubbles in immersion oil  | Remove bubbles   |
|   | Condenser aperture is closed or open too much                                       | Open or close properly   |
|   | Condenser is positioned too low   | Position the condenser upward  |

### ELECTRICAL PROBLEMS

| <b>Problem</b>                       | <b>Cause</b>        | <b>Solution</b>             |
|--------------------------------------|---------------------|-----------------------------|
| Lamp does not light when switched on | No electrical power | Check power cord connection |
|                                      | Lamp bulb burnt out | Replace bulb                |
|                                      | Fuse blown out      | Replace fuse                |

**IMAGE PROBLEMS**

| <b>Problem</b>             | <b>Cause</b>  | <b>Solution</b>   |
|----------------------------|---|---|
| Image moves while focusing | Specimen rises from stage surface                     | Secure the specimen in the slide holder                           |
|                            | Revolving nosepiece is not in the click-stop position | Revolve the nosepiece to the click-stop position                  |
| Image tinged yellow        | Lamp intensity is too low                             | Adjust the light intensity by rotating the intensity control dial |
| Image is too bright        | Lamp intensity is too high                            | Adjust the light intensity by rotating the intensity control dial |
| Insufficient brightness    | Lamp intensity is too low                             | Adjust the light intensity by rotating the intensity control dial |
|                            | Aperture diaphragm closed too far                     | Open to the proper setting  |
|                            | Condenser position too low                            | Position the condenser upward                                     |

**MECHANICAL PROBLEMS**

| <b>Problem</b>   | <b>Cause</b>                           | <b>Solution</b>                                 |
|--|--|---|
| Image will not focus with high power objectives                                    | Slide upside down                      | Turn the slide over so the cover glass faces up |
|  | Cover glass is too thick               | Use a 0.17mm cover glass                        |
| High power objective contacts slide when changed from low power objective          | Slide upside down                      | Turn the slide over so the cover glass faces up |
|  | Cover glass is too thick               | Use a 0.17mm cover glass                        |
|  | Diopter adjustment is not set properly | Readjust the diopter settings                   |
| Slippage of focus when using the coarse focusing knob<br>Fine focus is ineffective | Tension adjustment is set too low      | Increase the tension on the focusing knobs      |
|  | Tension adjustment is set too high     | Loosen the tension on the focusing knobs        |